

# **KOELIS**

## **Computer-Assisted Urology**

Puncture assistance




# Domain

Urologist	
Diagnosis	Therapy
	Prostate
	Kidney
	Bladder

- ✓ 3D tracking
- ✓ Image fusion
- ✓ Augmented reality
- ✓ Surgical Planning
- ✓ Measurement
- ✓ Surgical report
- ✓ ...

## ***Security-Confidence-Quality***

More accuracy  
 Less morbidity  
 Less radiation  
 Shorter stay  
 More insurance



# Products

- **Software**

- User friendly
- 2D/3D imaging & rendering
- Multimodal image fusion
- Computerized Surgical Protocol



- **Instruments**

- User friendly
- Sterilizable




- **Navigation platform**

- Stereo IR camera
- Tactile screen
- Control pedal

# Percutaneous Therapy of the Kidney

## CONTEXT

- Puncture needle
- Prone position
- CT scan / Fluo Rx
- Ultrasound probe + guide
- Local / Global anaesthesia



## PROBLEM

- Guide the needle fast & accurately

## SOLUTION


- CT navigation / Fluo navigation
- Real-time intra-operative planning

# Instrument and image localisation

- Infrared tracking of external references
  - On the calibrated C-arm, the patient and the needle guide
  - Spatial localisation of the fluoroscopic image and the needle trajectory in *patient space*




C-arm, needle guide and patient references



infrared camera and visualization system





RÉFÉRENCES

F		M
C	1	2
P		X

1	2	3
4	5	6
7	8	9

Paramètres d'images

-	☉	+
-	☀	+

Mémoriser Outil 1

Définir Axe

Vidéo directe du C-Arm




# Computerized Surgical Protocol

1. Install patient and patient reference on table
2. Calibrate C-arm with cage
3. Acquire 2-3 Fluo Rx (apnea)
4. Calibrate needle guide (with or without US)
5. Navigate virtual needle on Rx images
6. Perform the puncture confidently


# CT Navigation

- Hybrid localisation:
  - Tracking of needle and patient references
  - Image based registration of patient reference, visible in CT scan, with model of the reference shape



Needle reference

Camera





# Computerized Surgical Protocol

1. Install patient and patient reference on CT
2. Acquire 1 CT scan (apnea)
3. Calibrate needle guide (with or without US)
4. Match CT images onto patient reference
5. Navigate oblique CT slices
6. Perform the puncture confidently

# Work in progress

- 3DRT ultrasound based tracking
  - no external references
  - precise soft tissue tracking
  - patient reference replaced by reference volume of the target organ
  - volume stream registered with reference volume
  - deformation estimation possible
- Challenges
  - robustness
  - computational burden

# Navigation and robotics

- Common challenge of navigation systems and robotics: **Precise target localization**
  - Localization systems could be used for **robotic servoing**
    - Localizer based tracking
    - Image stream based tracking
    - Hybrid tracking



5, avenue du Grand Sablon  
BIOPOLIS  
38700 La Tronche, FRANCE  
[info@koelis.com](mailto:info@koelis.com)

